
Worthington Protease Protein Research Products





Primary Enzyme Producer

New Worthington Proteases for Protein Sample Prep

Trust your proteomics discoveries to Worthington. This catalog details our well-respected enzyme line for proteins requiring enzymatic digestion into peptides. Our expanded line of proteases enzymes appropriate for proteomic applications includes Arg-C, Glu-C, Lys-C and *SequENZ*[®] Modified Trypsin.

Sharing Our Enzyme Expertise

Education is a top priority at Worthington, and we have been authoring technical manuals and guides for over 70 years. Plus Worthington is cited in thousands of respected scientific journals across the globe.

To Support Your Research, We Provide Online Resources To Include:

- Advanced tissue search feature in the Tissue Dissociation Guide section of our website
- Comprehensive citations listings at: Bioz.com for Worthington products to obtain article snippets with technique filtering capabilities
- Collagenase Sampling Program to pre-test a particular lot of collagenase you are planning to use in your experiment. This free service allows you to pre-sample several different lots of collagenase at a time and select the best of the group for your application.

We invite you to work with us hand-in-hand to enhance our proteomics research technical library by submitting protocols, citations and articles referencing Worthington enzymes that can be shared with your colleagues. For details on submissions, contact your local Worthington Account Manager or forward suggestions to: techservice@worthington-biochem.com.



From writing the definitive Enzyme Manual decades ago, to our just released Introduction to Enzymes Video, our field account managers are here to assist you with the best products and services for your research.

Phone: 800.445.9603 • 732.942.1660 • Fax: 800.368.3108 • 732.942.9270 Worthington-Biochem.com ISO9001 Certified

SequENZ [®] Trypsin, Modified	3
Trypsin, Sequencing, Native	4
Endoproteinase-Arg-C	5
Endoproteinase Glu-C	6
Endoproteinase Lys-C	7
Chymotrypsin	8
Carboxypeptidase Y	9
Protease Products, Specifications and Applications Table	10-11
Ordering Information	12-13
Customer and Technical Support	14-16
International Distributors	17-19

Cataloa

Trypsin

Name

Source: Bovine Pancreas

CAS Number: 9002-07-7 I.U.B.: 3.4.21.4

Trypsin is a pancreatic serine protease with substrate specificity based upon positively charged lysine and arginine side chains. It is derived from a 34 kDa inactive precursor zymogen, trypsinogen, after enzymatic removal of an N-terminal 6-amino acid leader sequence resulting in the 23.8 kDa trypsin molecule. The optimum pH is 8.0. Trypsin is inhibited by organophosphorus compounds such as diisopropylfluorophosphate and natural inhibitors from pancreas. Soybean, lima bean, and egg white are also sources of natural inhibitors. Trypsin cleaves amide and ester bonds of Arg and Lys. The Worthington Sequencing Grade Trypsin has been further purified to remove trace contaminating proteases and autolysis products which could interfere in trypsin digestion experiments, and exhibits a single band on SDS PAGE.

Uses: For tissue culture work, Worthington trypsin, Codes: TRL, TRLS, TRLVMF and TRTVMF have been used by many researchers. Product Codes: TRSEQZ, TRSEQII and TRTPCK are typically used for protein sequencing, mapping and structure studies. Worthington modified sequencing grade trypsin, Product Code: TRSEQZ, is subjected to extensive purification to remove contaminating proteases and tryptic autolysis by-products which could affect the specificity of the digestion process. Subsequently, the enzyme is chemically modified to minimize the autolysis process as well as increase the stability. The modified trypsin is processed further to remove residual autodegradation products. The specificity of the enzyme is routinely checked after the chemical modification.

Stability/Storage: Most grades of Worthington trypsin are stable for 2-3 years when stored at 2-8°C. Protect from moisture.

Unit Definition: TAME Unit: One Unit hydrolyzes 1 micromole of p-toluene-sulfonyl-L-arginine methyl ester (TAME) per minute at 25°C, pH 8.2, in the presence of 10 mM calcium. One TAME Unit = 19.2 USP/NF units = 57.5 BAEE units.

Technical Notes: The Virus and Mycoplasma Free trypsin (Code: TRTVMF) has been filtered through an 0.22 micron pore size membrane, lyophilized, subjected to gamma irradiation, and tested for virus and mycoplasma.

Worthington certifies that all lots of Trypsin products are subjected to a pH of less than 3.0 for greater than five (5) hours during processing.

SequENZ® Trypsin, Modified,

Sequencing Grade

Trypsin, treated with L-(tosylamido-2-phenyl) ethyl chloromethyl ketone to inhibit contaminating chymotryptic activity, chemically modified to promote stability and further purified to remove autolysis fragments, resulting in a highly stable trypsin product resistant to autolysis while retaining specificity. Store at -20°C PROTECT FROM MOISTURE. **REQUIRES SPECIAL SHIPPING: ICE PACK** ≥ 150 Units per mg protein (≥ 8.625) **BAEE/2875 USP/NF** units per mg protein) LS02120 LS02122 LS02123 LS02124

4 x 25 µg

4 x 100 µg

1 ma

Bulk

Phone: 800.445.9603 • 732.942.1660 • Fax: 800.368.3108 • 732.942.9270 Worthington-Biochem.com ISO9001 Certified

TRSEQZ

Name	Activity	Catalog Number	Package	Code
Trypsin (Continued)				
Trypsin, Modified, Sequencing Grade, Solution Trypsin, treated with L-(tosylamido-2-phenyl) ethyl chloromethyl ketone to inhibit contaminating chymotryptic activity, chemically modified to promote stability and further purified to remove autolysis	≥ 150 Units per mg protein (≥ 8,625 BAEE/2875 USP/NF units per ml)	LS02150 LS02152 LS02155	250 μg 1000 μg Bulk	TRSEQZS
fragments, resulting in a highly stable trypsin product resistant to autolysis while retaining specificity. Supplied as a 1mg/ml solution. Store at -20°C REQUIRES SPECIAL SHIPPING: ICE PACK	. ,			
Trypsin, Purified, Sequencing Grade II Bovine trypsin that has been treated with L-(tosylamido-2-phenyl) ethyl chloromethyl ketone (TPCK) to inhibit contaminating chymotryptic activity and extensively purified to remove autolysis products. Supplied as a lyophilized powder. Store at -20°C. PROTECT FROM MOISTURE. REQUIRES SPECIAL SHIPPING: ICE PACK	≥ 150 Units per mg protein (≥ 8,625 BAEE/2875 USP/NF units per mg protein)	LS02115 LS02117 LS02119 LS02118	4 x 25 μg 4 x 100 μg 1 mg Bulk	TRSEQII
Trypsin, TPCK Treated A chromatographically purified, diafiltered, lyophilized powder that has been treated with L-(tosylamido-2-phenyl) ethyl chloromethyl ketone (TPCK) to inhibit contaminating chymotryptic activity (Kostka and Carpenter, <i>J. Biol. Chem.</i> <i>239</i> , 1799, 1964. Store at 2-8°C.	≥ 180 Units per mg protein (≥ 10,350 BAEE/3,450 USP/NF u/mg protein)	LS003740 LS003741 LS003744 LS003742	100 mg 500 mg 1 gm Bulk	TRTPCK

PROTECT FROM MOISTURE.

Clostripain (Endoproteinase-Arg-C)

Source: Clostridium histolyticum

I.U.B.: 3.4.22.8 CAS Number: 9028-00-6

Clostripain (Endoproteinase-Arg-C) is a two chain cysteine proteinase associated with collagenase and isolated from *Clostridium histolyticum*. It is highly specific for the carboxyl peptide bond of arginine. Clostripain is activated by dithiothreitol, cysteine, or other sulfhydryl containing reagents. The presence of calcium ions is essential. The enzyme is inhibited by oxidizing agents, divalent cations such as Co²⁺, Cu²⁺, Cd²⁺, and heavy metal ions. Citrate, borate, and Tris anions are less inhibitory.

Unit Definition: One Unit hydrolyzes one micromole of N-benzoyl-L-arginine ethyl ester per minute at 25°C, pH 7.6, in the presence of dithiothreitol.

		EQ
		СР
dry weight LS00164 LS00164	43 5 x 1 mg 46 10 mg	
)	Units per LS0016 Iry weight LS0016 LS0016	i protein LS02139 Bulk Units per LS001641 1 mg Iry weight LS001643 5 x 1 mg



Our mission is to provide superior tools from discovery research through larger scale bioprocessing applications.

ame	Activity	Catalog Number	Package	Code
rotease, Staphylococcus aureus (Er Source: Staphylococcus aureus V8	ndoproteinas	e Glu-C)		
I.U.B.: 3.4.21.19 CAS Number: 66676	6-43-5			
 side of either aspartic or glutamic acids. In the sites. It has a molecular weight of 27 kDa an <i>Staphylococcus aureus</i> V8 is inhibited by diis and NO₃⁻. Enzyme activity is determined by the 469 (1976). Stability/Storage: Autolysis occurs at temp Stable for 12 months at 2-8°C. 	d optimum pH of opropylfluorophos he casein digestio peratures greater t	4.0 and 7.8 with h phate and monov n assay describe han 40°C. The e	emoglobin as the substralent anions such as F- d by Drapeau, <i>Methods</i>	rate. Protease , Cl-, CH ₃ COO ⁻ <i>Enzymol., 45</i> , JSP 0.2% SDS.
Unit Definition: One unit causes a change of	51 0.001 A ₂₈₀ nm p	ber minute at 37 t	, pH 7.8 using casein a	s the substrate.
Protease, <i>S. aureus</i> Sequencing Grade Chromatographically purified according to Drapeau, <i>et al., J. Biol. Chem.,</i> 247, 6720 (1972). Supplied in vials containing 10 µg or 50 µg lyophilized powder for protein sequencing applications. Store at 2-8°C.	≥ 500 units per mg dry weight	LS02126 LS02128 LS02129	5x10 ug 5x50 ug Bulk	STSEG
Protease, S. aureus (Endoproteinase Glu	-			STA
Chromatographically purified according to	≥ 500 units	LS003608	1 mg	

per mg dry

weight

Drapeau, G., Boily, Y., and Houmard, J.,

J. Biol. Chem., 247, 6720 (1972).

A lyophilized powder. Store at 2-8°C.

Catalog

LS003605

LS003606

5 mg

Bulk



Worthington employees take pride in developing the highest quality enzymes.

Phone: 800.445.9603 • 732.942.1660 • Fax: 800.368.3108 • 732.942.9270 Worthington-Biochem.com ISO9001 Certified

LYSCSEQ

Endoproteinase Lys-C

Source: Lysobacter enzymogenes


I.U.B.: 3.4.21.50 CAS Number: 72561-05-8

Endoproteinase-Lys-C is a serine endoproteinase that specifically cleaves peptide bonds at the carboxyl side of lysine. It has a molecular weight of 30,000 daltons and is used in the optimal pH range of 7.0-9.0. Lys-C is inhibited by diisopropyl-fluorophosphate, TLCK, Aprotinin and Leupeptin.

Stability/Storage: Stable for 12 months at 2-8°C.

Unit Definition: One Unit will hydrolyze 1.0 micromole of N-p-tosyl-Gly-Pro-Lys p-nitroanilide per minute at pH 7.7 and 25C.

Endoproteinase Lys-C, Sequencing Grade			
Chromatographically purified.	≥ 150 units	LS02143	20 ug
A lyophilized powder supplied in	per mg	LS02144	5x20 ug
20 ug high recovery vials.	protein	LS02145	Bulk
Store at -20°C.			



Bringing you new ways to stay abreast of current methods, applications and products.

Chymotrypsin

Source: Bovine Pancreas

I.U.B.: 3.4.21.1 CAS Number: 9004-07-3

Chymotrypsin preferentially catalyzes the hydrolysis of peptide bonds involving L-isomers of tyrosine, phenylalanine and tryptophan. It also readily acts upon amides and esters of susceptible amino acids. Chymotrypsin catalyzes the hydrolysis of bonds of leucyl, methionyl, asparaginyl and glutamyl residues.

Stability/Storage:: The enzyme is stable for days in solution at pH 3.0 and for years as a dry powder at 2-8°C. Protect from moisture.

Unit Definition: One Unit hydrolyzes one micromole of benzoyl-L-tyrosine ethyl ester per minute at 25°C, pH 7.8 in the presence of calcium. An activity of 45 Units per mg using the above definition, is the equivalent of 10,000 optical density or 1330 N.F. units per mg using ATEE as a substrate.

1 BTEE unit = 29.5 USP/NF units.

Chymotrypsin, Alpha, TLCK Treated, Sequencing Grade Three times crystallized and treated with 1-chloro-3-tosylamido-7-amino-2-heptanone (TLCK) to inhibit trypsin activity (Shaw, <i>et al.,</i> <i>Biochemistry, 4</i> , 2219, 1965). Dialyzed against 1 mM HCl to remove autolysis products and low molecular weight contaminants. Supplied lyophilized in 25 ug and 100 ug high-recovery vials. Store at 2-8°C.	≥ 45 Units per mg protein	LS02130 LS02132	4 x 25 ug 4 x 100 ug	CDSEQ
Chymotrypsin, Alpha, TLCK Treated Three times crystallized and treated with 1-chloro-3-tosylamido-7-amino-2-heptanone (TLCK) to inhibit trypsin activity (Shaw, <i>et al.,</i> <i>Biochemistry, 4</i> , 2219 1965). Dialyzed against 1 mM HCI to remove autolysis products and low molecular weight contaminants. Supplied as a dialyzed, lyophilized powder. Store at 2-8°C.	≥ 45 Units per mg protein	LS001430 LS001432 LS001434 LS001438	25 mg 100 mg 1 gm Bulk	CDTLCK
Chymotrypsin, Alpha, Purified Chromatographically prepared by the procedure of Yapel <i>et al., J. Amer. Chem.</i> <i>Soc., 88,</i> 2573 (1966). A lyophilized powder. Store at 2-8°C.	≥ 45 Units per mg protein	LS001475 LS001479 LS001477	100 mg 1 gm Bulk	CDS
Chymotrypsin, Alpha, 3X Three times crystallized alpha chymotrypsin, which is an activation product of a three times crystallized zymogen. Dialyzed against 1 mM HCl and lyophilized. Store at 2-8°C.	≥ 45 Units per mg protein	LS001448 LS001450 LS001451 LS001453	250 mg 1 gm 10 gm Bulk	CDI
Chymotrypsin, Alpha, Crystallized Crystallized as zymogen and activated. Dialyzed against 1 mM HCl and lyophilized. Store at 2-8°C.	≥ 35 Units per mg protein	LS001333 LS001334 LS001332	1 gm 10 gm Bulk	CDAG

COY

Carboxypeptidase Y

Source: Yeast

I.U.B.: 3.4.16.5 CAS Number: 9046-67-7

Carboxypeptidase Y is a glycoprotein exopeptidase of the 'acid' and 'serine' class. Carboxypeptidase Y has a broad amino acid specificity, including proline and amidated amino acid residues.

Unit Definition: One Unit hydrolyzes 1 micromole of benzyl-oxycarbonyl-L-phenylalanyl-L-leucine per minute at 25°C, pH 6.5.

Carboxypeptidase Y

A highly purified preparation supplied	≥ 50 Units per	LS009070	1 mg
as a lyophilized powder.	mg protein	LS009068	5 mg
Store at -20°C.		LS009071	Bulk



From research and development to manufacturing, continuous quality improvement is everyone's job.

Enzyme	Specificity	Molecular Weight KDa	pH Ontimum	Extinction Coefficient	Common Substrates	Activators	Inhibitors	Product Code/ Applications
Drotocco Ea	· Protein Seque		Optimum	E1%, 280nm				
Carboxy- peptidase Y	H2-N-Rn-Y-†-X- COOH X ,Y= non-specific, prefers aromatic	64.0	4.5-6.0	15.0 (Hayashi <i>et al.</i> 1973, and Kuhn <i>et al.</i> 1973)	ATEE Bz-Phe-Ala-Leu Z-Phe-Ala	None required	APCK, Aprotinin DFP 4-Hydroxymercu- ribenzoate PMSF	COY, C-terminal sequencing & Modification/labeling of peptides and proteins
Chymotrypsin TLCK treated	-X-†-Y- X = aromatic Y = nonspecific	25.6	7.8-8.0	20.57 (Theoretical)	ATEE BTEE	None required	α-antitrypsin Aprotinin DFP, PMSF, TPCK α2-macroglobulin	CDSEQ, CDTLCK, Sequence analysis Peptide synthesis, mapping/finger- printing
Endo-Arg-C (Clostripain)	-Arg-†-Y- Y = nonspecific	53	7.4-7.8	16.57 (Theoretical)	BAEE	Ca ²⁺ Reducing agents	EDTA, TLCK, Tris Hg ²⁺ & other heavy metal ions	CPSEQ, CP, Peptide mapping & synthesis Sequence analysis Hydrolysis/ condensation of amide bonds
Endo-Glu-C (Staph. Protease V8)	-Glu-†-Y- (NH4 buffers pH 4, 7.8) -Asp-†-Y- (PO4 buffer pH 7.8)	27.0	4.0 & 7.8	4.26 (Houmard 1976)	Casein Z-Phe-Leu-Glu- 4NA	None required	DFP F-, CI-, Br-, CH3COO- NO3- α2-macroglobulin	STSEQ, STAP, Peptide mapping & sequence analysis
Endo-Lys-C	-Lys-+-Y- Y= nonspecific	30.0	7.0-9.0	18.63 (Theoretical)	N-p-Tosyl-Gly- Pro-Lys pNA	None required	DFP, TLCK, Aprotinin, Leupeptin	LYS-C, LYSEQ Peptide mapping and sequence analysis
SequENZ® Trypsin, Sequencing Grade, Modified								TRSEQZ, Modified Sequencing Grade, chemically modified to reduce autolysis Peptide mapping & sequence analysis
Trypsin, Sequencing Grade, Native	-X-†-Y- X = Arg, Lys Y = nonspecific	23.8	7.5-8.5	14.3	BAEE Casein TAME	Ca ²⁺ Lanthanide	Aprotinin, Benzamidine DFP, EDTA, Leupeptin α2-macroglobulin PMSF, TLCK Trypsin Inhibitors (egg white, lima	TRSEQII, Sequencing Grade, Native, Peptide mapping & sequence analysis
Trypsin, TPCK Treated							bean, pancreatic, soybean)	TRTPCK, TPCK Treated, Peptide mapping & sequence analysis Cleavage fusion proteins

L

Wort	hington Pr	otease	Produ	cts, Spec	ificatio	ns and A	Application	ns Table
Enzyme	Specificity	Molecular Weight KDa	pH Optimum	Extinction Coefficient E1%, 280nm	Common Substrates	Activators	Inhibitors	Product Code/ Applications
Partially Pur	ified for Tissue Di	issociation	and Proteir	n Digestion				
Collagenase	-Pro-X-†-Gly-Pro-Y- X = neutral Y = nonspecific	68-130	6.3-7.5	13.20 (CoIH, Theoretical) 13.40 (CoIG, Theoretical)	Collagen FALGPA Wünsch	Ca ²⁺ , Zn ²⁺	α 2-macroglobulin Cysteine, histidine DTT, 2-mercapto EDTA, EGTA Hg ²⁺ & other heavy metal ions <i>o</i> -phenanthroline	Collagenase products Tissue dissociation/ Primary cell isolation applications (see Tissue Dissocia- tion Guide for specific references)
Elastase	Elastin, -X-†-Y- X = uncharged, nonaromatic Y = nonspecific	25.9	8.0-8.5	21.8 (Theoretical)	Casein Denatured collagen Elastin, Fibrin Suc-Ala3-NA	None required	α -antitrypsin DFP α 2-macroglobulin PMSF	ES/ESL, suspension/ lyo powder, Tissue Dissociation/ Primary cell isolation applications (see Tissue Dissocia- tion Guide for specific references)
Neutral Protease (Dispase®)	-X-†-Leu/Phe-†-Y- X/Y = nonspecific	36.0	5.9-7.0	13.96 (Theoretical)	BAEE Casein	Ca^{2+} , Mg^{2+} , Mn^{2+} , Fe^{2+} , and Al^{3+}	EDTA, EGTA Hg ²⁺ & other heavy metal ions <i>o</i> -phenanthroline	NPRO/NPR02, Tissue Dissociation/ Primary cell isolation and cell harvesting applications (see Tissue Dissociation Guide for specific references)
Papain	-X-†-Y- X = nonspecific but Arg, Lys and Phe preferred Y = nonspecific	23.0	6.0-7.0	22.88 (Theoretical)	BAEE	Cysteine EDTA Reducing agents GSH, NBS	AEBSF, Antipain Cystatin, Leupeptin α 2-macroglobulin Hg ²⁺ & other heavy metal ions DFP, PMSF TLCK, TPCK, E-64	PAP/PAPL, suspen- sion/lyo powder, Neural tissue dissociation/ primary cell isolation applications (see Tissue Dissociation Guide for specific references) Antibody cleavage RBC modification
Pepsin	$-X-\uparrow-Y-$ X = nonspecific but aromatic & hydro- phobic preferred Y \neq Ala, Gly, Val	34.6	1.0-4.0 unstable ≥5	14.39 (Theoretical)	Casein Hemoglobin	None required	Pepstatin A Diazoketones Epoxides	PM, Collagen bioprocessing/ purification Antibody cleavage
Proteinase K	-X-†-Y- X = nonspecific but aliphatic, aromatic & hydrophobic preferred Y = nonspecific	28.9	7.5-12	12.6 (Theoretical)	Casein Hemoglobin Keratin	Ca ²⁺ Active in 0.5- 1% SDS	DFP EGTA PMSF	PROK, PROKS, DNA/RNA purification
Trypsin	-X-†-Y- X = Arg, Lys Y = nonspecific	23.8	7.5-8.5	14.3	BAEE Casein TAME	Ca ²⁺ Lanthanide	Aprotinin, Benzamidine DFP, EDTA, Leupeptin α2-macroglobulin PMSF, TLCK Trypsin Inhibitors (LBI, OI, SI/SIC)	Trypsin products Protein Digestion/ Sequencing (purified) Tissue dissociation/ Primary cell isolation applications (see Tissue Dissociation Guide for specific references)

Worthington provides various options to make ordering fast and convenient:

• Call:	1.800.445.9603 (8am – 5:30pm EST Mon. – Fri.) 1.732.942.1660
• Fax:	1.800.368.3108 1.732.942.9270 (24hr/day, 7 days/week)
• E-mail:	custservice@Worthington-Biochem.com techservice@Worthington-Biochem.com
Online/Website:	Worthington-Biochem.com TissueDissociation.com
• Write:	Worthington Biochemical Corporation 730 Vassar Avenue Lakewood, New Jersey 08701 U.S.A

Terms of Sale

Prices effective January, 2021 and are subject to change without notice. Additional charges added for shipping, which will be prepaid and added to the invoice unless other arrangements are made at the time of ordering. Insurance will be charged for higher-value shipments at our discretion.

Payment terms are Net 30 Days, F.O.B. Origin, Lakewood, New Jersey USA, payable in US dollars. All checks must be drawn on a US bank or payment made by wire transfer. Past due accounts may be charged a 1.5% per month late payment fee.

Complete Standard Terms and Conditions of Sale available on our website.

VISA, MasterCard & American Express are also accepted.



Quantity	Price
1 to 4	List Price
5 to 9	5% off list
10 to 19	10% off list
20 or more	15% off list
Bulk	Inquire

Quantity Discounts

Standing Orders & Additional Discounts

For orders of greater than 25 packages, or orders of material packed in bulk, contact your representative or the Bulk Sales Office for special pricing consideration. Standing orders may also qualify for discounts. We welcome long-term use projections for which we can consider special rates. Large institutional buyers should contact their representative regarding special purchasing agreements.

Bulk, Contract/Custom & OEM Enzyme Purchasing

As a primary manufacturer, Worthington can supply products in a wide range of purity and activity specifications and in large-scale bulk quantities at substantial discounts. In addition, we welcome inquiries for contract and custom manufacturing, custom analysis, and special packaging for OEM applications. Several products are listed as Bulk Only in this catalog due to limited availability. Please contact Customer Service or our Bulk Sales Office to discuss your specific requirements.

Complete Standard Terms and Conditions of Sale available on our website.

Available 8:00 AM to 5:00 PM Eastern Time Monday through Friday. We can be contacted 24 hours a day by fax, e-mail or through our website.

Worthington makes the products we sell and welcome your questions and suggestions. Because we are a primary manufacturer we have ready access to all production and quality control records of our products by lot number.

Our years of experience in enzyme purification put us in a position to assist individual researchers with special needs. We frequently do customized preparations of entirely new products. We can make modifications of a regular production procedure on a custom basis. Furthermore, our quality control department can do special testing if needed.

Online Sampling Program



Need help with protocols? Ask a representative how we can help update you with our latest technical tools.

Our position as the principal manufacturer of research grade collagenase makes possible our Collagenase Sampling Program. Under the program, we provide 100 mg samples of up to three different lots of collagenase for evaluation in your own cell isolation systems. A period of 60 days is allowed for your evaluation of these samples. A minimum of 3 grams of each lot of collagenase will be placed on HOLD, reserved in your name. When you determine

which lot performs best for you, specify the lot desired when ordering. The only requirement, once a suitable lot of collagenase is found, is that you purchase a minimum of 3 grams of the material.

There is no charge for participating in the Collagenase Sampling Program. Contact your representative or our Technical Service group for more information: techservice@Worthington-Biochem.com.

Collagenase Lot Selection Tool

Worthington's Collagenase Lot Selection Tool is available online at our website. This feature was designed to help researchers select and evaluate current collagenase lots that match previous lots or desired activity profiles. Users may enter target values for collagenase, caseinase, clostripan, and tryptic activities or specify previous lot numbers. Each value can be weighted based on the relative level of importance to the application. After the search for matches is completed, a ranked list of collagenase lots currently available is generated. The selected lots can then be sampled simply by using the built-in link to the Free Collagenase Sampling Program. As always, Worthington Customer and Technical Service personnel are available via phone and e-mail to assist with collagenase or any other products.

Worthington Biochemical Corporation is company-wide ISO9001 certified and operates according to GMP guidelines. Our initial ISO assessment audit was performed by ANAB-accredited SGS US Testing Company, Systems & Services Certification in 2005 with continuous successful re-certifications.

Product Use

All Worthington products are sold for manufacturing, research, and laboratory use only by properly trained and authorized personnel. Researchers and clinical laboratory personnel intending to use any of these products for medical investigation on humans are solely responsible for such use, and for compliance with the pertinent regulations of the United States Food & Drug Administration (USFDA) and other regulations. We do not assume liability for damages resulting from the use of these products or from their use in violation of patent or other rights.

U.S.D.A. Certified Raw Materials

All products from animal sources are produced from starting materials of United States Department of Agriculture (USDA) or equivalent approved origin, collected in USDA or equivalent approved facilities, inspected to be free of disease and suitable for exportation. Certificates of Origin are available upon request.

Animal Free Products

Several Animal Free (AF) nucleases, proteases and other products are also available to eliminate BSE/TSE and mammalian viral risks. Please inquire. All animal free products are designated with this symbol for ease of use.

Product Returns

Authorization for any product return must be obtained from Worthington Biochemical Corporation (Customer Service Department), or its authorized representative, prior to the return of product. This authorization is required to insure the proper return of material and, if applicable, the correct issuance of credit. There is no provision for credit of misused, improperly stored or outdated material. Product(s) must be returned in the same condition as received within 30 days of the original shipment by Worthington Biochemical Corporation. A restocking fee may be charged.

Complete Standard Terms and Conditions of Sale available on our website.

Hundreds of pages of technical, product and reference information are available online at: Worthington-Biochem.com

Worthington Enzyme Manual Tissue Dissociation Guide Catalog and Price List Online Collagenase Lot Selection Tool

Additional Features

Complete Searchability Current Collagenase Lot Availability and Activities Online Ordering and Technical Service Updated Announcements and Exhibit Schedules International Distributor Listing Useful Links

Linked in

Visit us at: facebook.



Addressing your individual needs, we value every customer interaction. Let us know how we are doing.

Phone: 800.445.9603 • 732.942.1660 • Fax: 800.368.3108 • 732.942.9270 Worthington-Biochem.com ISO9001 Certified

Name

Australia:

ScimaR Email: scimar@bigpond.net.au

Austria:

Eubio

Email: eubio@eubio.at Website: www.eubio.at

Belgium:

GESTIMED s.p.r.l/b.v.b.a.

Email: info@gestimed.be Website: www.gestimed.be

Brazil:

Sellex Email: vendas@sellex.com Website: www.sellex.com

Sinapse Biotecnologia Ltda

Website: www.sinapsebiotechnologia.com.br

Canada:

Cedarlane

Email: info@cedarlanelabs.com Website: www.cedarlanelabs.com

Chile:

Fermelo S.A. Email: contacto@fermelo.cl Website: www.fermelo.cl

China:

Activity

4A Biotech Company, Ltd.

Package

Email: info@4aBio.com Website: www.4aBio.com

Cataloa

Number

EQUL Ltd.

Email: info@equl.com Website: www.equl.com

Gene Company Ltd.

Email: info@genehk.com Website: www.genehk.com

Genetimes Excell Technology, Inc.

Email: order@genetimes.com.cn Website: www.genetimes.com.cn

MACGENE (Beijing) Biotechnology Ltd.

Email: order@macgene.com Website: www.macgene.com

Shanghai 78 Biosciences Co., Ltd.

Website: www.qfbio.com

Shanghai Universal BioTech Co., Ltd.

Email: purchasing@univ-bio.com Website: www.univ-bio.com

Denmark:

BioNordika (Denmark) A/S

Email: info@bionordika.dk Website: www.bionordika.dk

England/UK:

Lorne Laboratories, Ltd. Email: info@lornelabs.com Website: www.lornelabs.com

Finland:

BioNordika Finland OY

Email: info@bionordika.fi Website: www.bionordika.fi

Sigma-Aldrich Finland OY

Email: NordicOrder@sial.com

Activity

France:

Name

Serlabo Technologies

Email: baudet@serlabo.fr Website: www.serlabo.fr

Germany:

CellSystems Biotechnologie Vertrieb GmbH Email: info@cellsystems.de

Website: www.cellsystems.de

Hong Kong:

Gene Company, Ltd.

Email: info@genehk.com Website: www.genehk.com

Line Analytics Life Sciences, Ltd.

Email: sales.lals@hcdh.com.hk Website: www.hcdl.hcdh.com.hk

India:

Arun & Company

Email: diagnostics@zytex.com Website: www.arunandco.com

Rahesh & Company

Email: diagnostics@zytex.com Website: www.arunandco.com

Indonesia:

CV Gamma Scientific Biolab

Email: gammascientific@gmail.com Website: www.gammascientificbiolab.co.id

PT Genetika Science Indonesia

Email: info@ptgenetika.com Website: www.ptgenetika.com

Israel:

Enco Scientific Services, Ltd.

Email: info@enco.co.il Website: www.enco.co.il

Italy:

D.B.A. Italia, s.r.l. Email: info@dbaitalia.it Website: www.dbaitalia.it

Japan:

Funakoshi Company, Ltd. Email: reagent@funakoshi.co.jp

Website: www.funakoshi.co.jp

Korea:

Chayon Laboratories, Inc. Email: info@chayon.co.kr Website: www.chayon.co.kr

Dong In Biotech Co., Ltd.

Email: info@donginbio.com Website: www.donginbio.com

Kim & Friends, Inc.

Email: kimnfriends@hanmail.net Website: www.kimnfriends.co.kr

LRS Laboratories, Inc.

Email: info@lrslab.co.kr Website: www.lrslab.co.kr

Malaysia:

BioSynTech Malaysia Group Sdn Bhd Email: sales@bstmgroup.com Website: www.bstmgroup.com

Essen-Haus Sdn. Bhd.

Email: sales@essen-haus.com.my Website: www.essen-haus.com.my

i-DNA Biotechnology (M) Dsn Bhd

Email: sales@i-dna.com.my Website: www.i-dna.com.my

Mexico:

Netherlands:

Interchema C.V. - Antonides Email: info@interchema.com Website: www.interchema.com

Norway:

Nerliens Meszansky AS Email: info@nmas.no Website: www.nmas.no

Singapore:

i-DNA Biotechnology Pte Ltd Email: info@i-dna.sg Website: www.i-dna.sg

Spain:

LabClinics, S.A. Email: info@labclinics.com

Sweden:

BioNordika (Sweden) AB Email: info@bionordika.se Website: www.bionordika.se

Switzerland:

BioConcept

E-mail: info@bioconcept.ch Website: www.bioconcept.ch

Taiwan:

Integrated Bio LTD Email: info@integrated-bio.com Website: www.integrated-bio.com

Johnson & Annie Company, Ltd.

Email: home@johnsonandannie.com Website: www.johnsonandannie.com

Protech Technology Enterprise Co. Ltd.

Email: service@bio-protech.com.tw Website: www.bio-protech.com.tw

Thailand:

Pacific Science Co., Ltd.

Email: pacscien@ksc.th.com Website: www.pacificscience.co.th

Please visit: Worthington-Biochem.com For territories not listed or contact International Sales, Email: international@Worthington-Biochem.com

Our Driving Force is Quality

The energy that drives all Worthington employees is the deep commitment to consistently produce superior quality products while providing personal service.

Our ISO9001 certified enzyme manufacturing facility and proprietary processes are designed specifically for the life science community, bioprocessing and OEM/private label business partners. We are proud to be a primary producer with over 90% of our products manufactured at our core facility located in Lakewood, New Jersey, USA.

Value Verified

 \checkmark

ISO9001 Certified Quality Management System All Processes Fully Documented and Traceable from Raw Materials through Final Shipment

 \checkmark

Biopharmaceutical and Bioprocessing Respected Manufacturer of High Quality, Lot-to-Lot Consistent Enzymes for Over 70 Years

Primary Producer For Biotech, Life Science Research, Diagnostic,

 \checkmark

Multi-Scale Extraction, Fermentation, Protein Purification, Lyophilization and Packaging Capabilities Produced Under GMP Guidelines

Internal Testing Capabilities include Enzyme, Protein, and Related Biochemical Analysis and Characterization

 \checkmark

Animal-Free Certified and USDA Approved Animal-Sourced Products Suitable for Worldwide Exportation

 \checkmark

Flexible and Responsive Production Scheduling Sample Lots and Bulk Inventories Available for Immediate Shipment



Superior Customer Services and Technical Support — Before, During and After Sale

Volume Purchases, Bulk Packaging, Standing Orders, OEM/Supply Agreements, and Special Arrangements



730 Vassar Avenue • Lakewood, New Jersey 08701



Worthington-^ • Tel: 800.445.9603 • 732.942.1660 • Fax: 800.368.3108 • 732.942.9270