	Subscribe (/subscribe)	About (/about)		Welcome,	, Jim. You have trial access to		Learn More	
	Editorial Boards (/editorial	-boards) Auth	nors (/publish)				<u>(/access)</u> .	
	Blog (http://blog.jove.com)							
<u>(/)</u>		Separate s	search terms	with commas	Search	Browse -	Sign Out	

# **Primary Culture of Adult Rat Heart Myocytes**

Xianghua Xu (/author/Xianghua\_Xu)<sup>1</sup>, Henry M. Colecraft (/author/Henry+M.\_Colecraft)<sup>1,2</sup>

<sup>1</sup>Department of Physiology and Cellular Biophysics, Columbia University (/institutions/NA-north-america/US-united-states/NY-new-york/24811-columbia-university), <sup>2</sup>Department of Pharmacology, Columbia University (/institutions/NA-north-america/US-united-states/NY-new-york/24811-columbia-university)



Video Article Chapters	
0:04	Title
0:52	Introduction
1:18	Preparation for Surgery
4:06	The Rat Heart Dissection
5:05	Myocytes Isolation
10:15	Myocyte Cell Culture
10:56	Representative Outcome: Cultured Myocytes
11:20	Conclusion



#### **Related Videos**



Isolation and Genetic Manipulation of Adult... Published 9/17/2009

Isolation of Rat Portal

Fibroblasts by In... Published 6/29/2012



Protocol for Culturing Sympathetic Neurons from... Published 1/30/2009

Exploring Arterial

Smooth Muscle Kv7 Potassium... Published 9/14/2012

#### Download

EndNote, RefWorks, Reference Manager, ProCite (/legacy/GetRIS.php?id=1308)

JoVE PDF Article (/pdf/1308)

Materials List (/pdf-materials/1308)

XML File (/pdf/default.aspx? ID=1308)

#### Information

Date Published	6/16/2009
JoVE Section	General (/general)
JoVE Issue	28 (/archive/28)
Comments	6 Comments
View Count	35,125 (Details)
DOI	doi:10.3791/1308

#### **Cite this Article**

Xu, X., Colecraft, H. M. Primary Culture of Adult Rat Heart Myocytes. J. Vis. Exp. (28), e1308, doi:10.3791/1308 (2009).

# Abstract

Cultured primary adult rodent heart cells are an important model system for cardiovascular research. Nevertheless, establishment of robust, viable cultured adult myocytes can be a technically challenging, ratelimiting step for many researchers. Here we described a protocol to obtain a high yield of adult rat heart myocytes that remain viable in culture for several days. The heart is isolated and perfused with collagenase and protease under low Ca<sup>2+</sup> conditions to recover single myocytes. Ca<sup>2+</sup>-tolerant cells are obtained by stepwise increases in extracellular  $Ca^{2+}$  concentration in three subsequent wash steps. Cells are filtered, resuspended in culture medium, and plated on laminin coated slips. Cultured myocytes obtained using this protocol are viable for up to four days and are suitable for most experiments including electrophysiology, biochemistry, imaging and molecular biology.

Protocol

#### Keywords

#### Part I. Preparation for surgery

- The day before the rat surgery, make certain that all solutions are made, but do not add any enzymes yet. Buffers A, B and the wash buffer should all be filtered for sterility prior to storage. The buffers should all be stored at 4 °C.
- 2. 6-well tissue culture plates must be plated with 10-20  $\mu$ g/mL laminin in 2 mL 1XSFM the day before as well, and should be stored at 37 °C in a CO<sub>2</sub> incubator.
- 3. On the day of the surgery, prepare the surgical instruments by sterilizing the clamp, scissors, and forceps with 70 % ethanol and let them dry on paper towel.
- 4. A 1 mL syringe filled with 0.5 mL Heparin Solution (90 U/ml) is made ready.
- 5. Now it is good practice to wash the perfusion apparatus by running through 70% ethanol, using the peristaltic pump in reverse. Once finished rinsing with ethanol, rinse the apparatus with double distilled water, and finally rinse with Buffer A. The flow through is collected in a beaker placed in a 37 °C water bath.
- 6. To prepare the cleaned perfusion apparatus, fill the right syringe tube with 40 mL of Buffer A and the left syringe tube with 40 mL of Buffer B.
- Prime the perfusion tubing to the tip of the catheter by allowing Buffer A to flow through the apparatus.
  Refill to the buffer A syringe to the 40 mL level. Make sure there are no air bubbles trapped in the tubing after this step.
- 8. Now, close the buffer A syringe stopcock valve and repeat the process with Buffer B such that the perfusion tubing is now primed with Buffer B.
- 9. Finally, the perfusion apparatus is left with the Buffer B stopcock valve open, but with the flow stopped using the regulator on the IV line.
- 10. Discard the buffer flow-though in the collecting beaker.
- 11. The last thing to do before the surgery is to add the required enzymes to your stock solution, and then filter the enzyme solution just as the other solutions were filtered the day before. Once filtered, this solution can be left at room temperature.

#### Part II. The Rat heart dissection

- 1. Following standard protocol euthanize the rat with Isoflurane in a gas chamber.
- 2. Sterilize the incision area on the thorax with 70% ethanol.
- 3. Open chest using scissors and expose heart.
- 4. Inject left ventricle with 0.5 mL Heparin solution (90 U/mL).
- 5. Remove heart with large portion of aorta intact, place in ice cold Buffer B.

#### Part III. Myocytes isolation

- A typical rat heart should yield a high fraction of rod-shaped myocytes (as opposed to dead rounded cells) if the following procedure is properly adhered to.
- 2. To begin, clamp aorta to the catheter. The tip of the catheter should not be pushed too far into the heart to ensure good perfusion of the heart through the coronary artery.
- 3. Once clamped, tie the aorta to the catheter with a suture.
- 4. Next, open the IV line regulator to allow a fast drip rate (~60 drops/min) and allow Buffer B to wash out

Cellular Biology (/keyword/cellular+biology), Issue 28 (/archive/28), Rat (/keyword/rat), heart (/keyword/heart), isolation (/keyword/isolation), primary culture (/keyword/primary+culture), cardiology (/keyword/cardiology)

Translate this page to:

Choose Language...

#### **Navigate this Article**

Cite this article Abstract Protocol Discussion Disclosures Materials References Ask the author Top of page the heart.

- 5. During Buffer B wash, use the small scissors and forceps to remove the atria and any fatty or lung tissue clinging to the heart.
- 6. After Buffer B is finished, switch the flow to Buffer A.
- 7. While Buffer A is allowed to flow for 5 minutes, several small tasks should be quickly performed.
- 8. First, warm the Enzyme Solution in a 37 °C water bath.
- Finally, when Buffer A is depleted from the syringe (but still present in the tubing), load 50 mL of the Enzyme Solution in syringe A of the perfusion apparatus
- 10. It is now necessary to discard all the flow-through from the collection beaker in the water bath (by pipetting), until the Enzyme Solution perfuses heart.
- 11. As soon as the Enzyme Solution perfuses the heart, activate the peristaltic pump which is set up to transfer the enzyme solution from the collecting beaker to replenish syringe A.
- 12. Allow the enzyme solution to flow through the heart for 10 min. As the heart digests, it begins to look bloated.
- 13. Add 37.5  $\mu$ L of 0.1 M CaCl<sub>2</sub> to the enzyme solution in syringe A to give an effective concentration of 0.1 mM Ca<sup>2+</sup>.
- 14. After 10 minutes increase the concentration of calcium to 0.2 mM by adding an additional 50  $\mu$ L of 0.1 M CaCl<sub>2</sub> to syringe A and let the perfusion proceed for a further 10 minutes.
- 15. Cut off the ventricles and transfer to a small sterile beaker containing 20 mL Enzyme Solution.
- 16. In the beaker increase the Calcium concentration to 0.4 mM by adding 40 more  $\mu L$  of 0.1 M CaCl\_2.
- 17. Gently mince the heart into 10 or more pieces with a pair of small sterile scissors.
- 18. Incubate the beaker for 5 minutes at 37 °C with gentle rocking.
- Add 40 mL of 0.1 M CaCl<sub>2</sub> to give effective concentration of 0.6 mM Ca<sup>2+</sup> and gently triturate the heart pieces with a plastic transfer pipet 3-5 times.
- 20. After incubating for an additional 5 minutes at 37 °C, add 40 mL of 0.1 M CaCl<sub>2</sub> and gently triturate as before.
- Separate digested single myocytes from non-digested connective tissue using filtration through a sterile 500 μm mesh.
- 22. Allow the cells to settle in a 50 ml tube for 10 minutes at room temperature.
- 23. Discard supernatant with a transfer pipet.
- 24. Gently resuspend cells in Wash Buffer #1 and allow the cells to settle for 10-20 min at room temperature.
- 25. While the cells are settling, take a small aliquot, and use this time as an opportunity to assess the quality and viability of the cells in suspension under an inverted microscope. A good preparation will have a high proportion (>80%) of rod like cells with crisp striations.
- 26. Once the cells have settled, discard the supernatant and gently resuspend cells in Wash Buffer #2. Let the cells settle at room temperature, which should again take 10-20 minutes, and discard the supernatant.
- 27. Repeat the washing process once more with Wash Buffer #3, and the cells will be ready for culture.

#### Part IV. Myocyte cell culture

- Gently resuspend cells in 5% Serum Media. Before transferring the cells to tissue culture plates, wash the plates with 1x SFM. Transfer the cells at a desired density and incubate in a CO<sub>2</sub> tissue culture incubator for 3-5 hours.
- 2. Switch the media to 1x SFM.

3. Cells can be cultured for up to four days and used as required for experiments.

# Part V. Representative Results/Outcome

There should be more than 70% live heart myocytes under inverted microscope when the protocol is done correctly.



Table 1: Solution recipe for 1 L 10X Ca<sup>2+</sup>-free KH Buffer:

	mass	Final Concentration
NaCl	68.96 g	1180 mM
КСІ	3.58 g	48 mM
HEPES	59.58 g	250 mM
K <sub>2</sub> HPO <sub>4</sub>	2.85 g	12.5 mM
MgSO <sub>4</sub>	3.08 g	12.5 mM

Table 2: Solution recipe for 1 L Buffer A:

Add 1.98 g glucose to 100 mL 10X KH buffer and bring to a final volume of 1 L. Adjust osmolarity to near physiology condition (~300) with glucose. Bring to pH 7.4 with NaOH.

-	Mass	Final Concentration			
10X KH Buffer	100 mL	100 ml			
Glucose	1.98 g	11 mM			
Bring to 1 L					
Bring to pH 7.4 with NaOH					
Adjust osmolarity to near physiological condition(~300) with glucose.					

### Table 3: Solution recipe for Enzyme Buffers

	mass	Final Concentration
Collagenase Type II	25 mg	0.05% (w/v)
Protease XIV	10 mg	0.02% (w/v)

BDM	0.025 g	5 mM
Carnitine	0.020 g	2 mM
Taurine	0.031 g	5 mM
Glutamic Acid	0.020 g	2 mM
0.1 M CaCl <sub>2</sub>	12.5 uL	25 uM
Buffer A		Fill to 50 ml

# Table 4: Solution recipe for 25X BDM/Taurine/BSA (B/T/B)

	mass	Final Concentration
BDM	0.076 g	125 mM
Taurine	0.094 g	125 mM
BSA		0.1% (w/v)
Buffer A		Fill to 6 mL

#### Table 5: Solution recipe for Buffer B and Wash Buffers

Buffer:	В	#1	#2	#3
Add:				
0.5 M CaCl <sub>2</sub>	0.4 mL	0.1 mL	0.125 mL	0.15 mL
Final concentration of CaCl <sub>2</sub>	1 mM	1 mM	1.25 mM	1.5 mM
25X B/T/B		2 mL	2 mL	2 mL
Buffer A to final volume	200 mL	50 mL	50 mL	50 mL

Table 6: Solution Recipe for Culturing Media

	Final Concentration
	(after dilution)
0.099 g	10 mM (5 mM)
0.063 g	10 mM (5 mM)
0.075 g	10 mM (5 mM)
0.5 ml	1% (0.5%) (v/v)
45 mL	Fill to 50 mL
	Final Concentration
	0.099 g 0.063 g 0.075 g 0.5 ml 45 mL

2X SFM	25 mL	1X
Media 199		Fill to 50 mL
1X 5% Serum Media		
		Final Concentration
2X SFM	25 mL	1X
2X SFM FBS	25 mL 2.5 mL	1X 5% (v/v)

### Discussion

A critical step is the speed with which the isolated heart is hung up on the perfusion system. The length of the enzymatic digestion period may be a little different for each rat. The adjustment depends on how soft the heart becomes after the regular period of digestion. The slowly recovery of  $Ca^{2+}$  after enzyme digestion is essential for obtaining  $Ca^{2+}$ -tolerant healthy cells.

For guinea pig, the protocol is the same except hyaluronidase is used instead of Protease XIV. Typically, we find the that initial percentage of living cells is lower for guinea pig compared to rat, although they survive just as long in culture.

#### Disclosures

#### Materials

Name	Company	Catalog Number	Comments
NaCl	Sigma-Aldrich	S6191	
КСІ	Sigma-Aldrich	P9541	
HEPES	Sigma-Aldrich	H4034	
K <sub>2</sub> HPO <sub>4</sub>	Sigma-Aldrich	P9666	
MgSO <sub>4</sub>	Fluka	63068	
Glucose	Sigma-Aldrich	G7021	
CaCl <sub>2</sub>	Sigma-Aldrich	C7902	
Heparin Sodium Salt	Sigma-Aldrich	H4784	
Collagenase Type II	Worthington Biochemical	LS004176	
Protease XIV	Sigma-Aldrich	P5147	
2,3-Butanedione Monoxime (BDM)	Sigma-Aldrich	B0753	
Carnitine	Sigma-Aldrich	C9500	
Taurine	Sigma-Aldrich	T0625	

Glutamic Acid	Sigma-Aldrich	G1501	
BSA	Sigma-Aldrich	A3294	
Creatine	Sigma-Aldrich	C3630	
Antibiotic	Cellgro	30-009-CI	
Media 199	GIBCO, by Life Technologies	11150	
FBS	Hyclone	SH30071.03	
Laminin	BD Biosciences	354232	
water bath	Fisher Scientific	15-462-5	
variable flow chemical pump	Fisher Scientific	15-077-67	
6-well culture plate	Corning	3516	

#### References

- 1. Miriyala J, Nguyen T, Yue DT, & Colecraft HM. Role of CaVbeta subunits, and lack of functional reserve, in protein kinase A modulation of cardiac CaV1.2 channels. Circ Res.;102(7):e54-64 (2008).
- 2. Takahashi SX, Mittman S, & Colecraft HM. Distinctive modulatory effects of five human auxiliary beta2 subunit splice variants on L-type calcium channel gating. Biophys J.;84(5):3007-21 (2003).

### Ask the Author

# **6** Comments

Can this protocol used for mouse hearts?

# 1 (/video/1308/primary-culture-of-adult-rat-heart-

It works fine for rat and guinea pig and we haven't had a chance to try it on mouse.

### 1.1 (/video/1308/primary-culture-of-adult-rat-heart-

I have a question. Can I use collagenase type I instead of type II ??

### 2 (/video/1308/primary-culture-of-adult-rat-heart-

I have a question. Can I use collagenase type I instead of type II ??

# 3 (/video/1308/primary-culture-of-adult-rat-heart-

Some works employ just left ventricle for myocyte isolation but here, we read "ventricles" through the text. Should I assume that both left and right ventricles ere used for isolation?

# 4 (/video/1308/primary-culture-of-adult-rat-heart-

Yes, both left and right ventricles would be present in the isolated myocytes.

### 4.1 (/video/1308/primary-culture-of-adult-rat-heart-

Thank you for your kind reply, Dr. Colecraft.

I have another two questions:

a) how many cells does this procedure yield?

b) how many cells do you plate, in a 6-well plate?

Thanks a lot!

#### 5 (/video/1308/primary-culture-of-adult-rat-heart-

can you please tell me why a few minutes after I start the heart reperfusion with Collagenase Type II, the heart to turn brownish.

Then when I incubate the heart at 37C with gentle rocking, I don't have a lot of myocytes after the filtration.

Thank you so much

### 6 (/video/1308/primary-culture-of-adult-rat-heart-

### Post a Question / Comment / Request

Your comment (in English, please):

Allowed HTML tags: i, b, u, sup, sub

You are signed in as jzacka@worthington-biochem.com

# **About JoVE**

The Journal of Visualized Experiments (JoVE) is a PubMed-indexed video journal. Our mission is to increase the productivity of scientific research.

**Follow JoVE** 



blog (http://blog.jove.com) rss (/rss)

Post Comment

twitter (http://twitter.com/JoVEjournal) newsletter

Subscribe (/subscribe)	About (/about)	Editorial Boards (/editoria	ll-boards)	Publish (/publish)	Press (/about/press-releases)
Advertise (http://industry.jove.com	) JoVE Shop	(http://shop.jove.com/)	Contact (/a	about/contact#contact)	Troubleshooting (/about/contact)

Copyright © JoVE 2006-2013. All Rights Reserved. | Policies (/policies) | License Agreement (/files/documents/JoVE\_library\_license\_agreement.pdf) | ISSN 1940-087X MyJoVE Corp. | 17 Sellers St. | Cambridge, MA | 02139 | tel. 617.945.9051 (tel:617.945.9051)